

Accumen Performance Suite 2.1.3



Product Release Notes

Accumen Performance Suite®
2.1.3
May 2023

Product Detail and Version

Product Name	Version	Date	Comment
Accumen Performance Suite	2.1.3	May 2023	

Revision History

Date	Version	Description	Author
05/11/2023	1	Release Notes	Celena Blackwood

Feature Release list contained in Version 2.1.3

New + Improved

1. APS Menu Enhancements
 - a. Breadcrumb
 - b. Click Path
2. Excel file Download
3. Laboratory Quality and Service Scorecard and Detail Display
4. Laboratory Quality and Service Detail Date Range
5. Laboratory Orders Query TAT Metric Update
 - a. Order to Collect TAT
6. Laboratory CPT Detail
 - a. Organization of display
 - b. Chart title during filtering
7. Message Center (Defect)
8. Laboratory Quality and Service Detail Filtering (Defect)

Detailed Product Release Description

1. APS Menu Enhancements

Addition of Breadcrumbs to Toolbar

Breadcrumb has been added to the top left of screen displaying the users click path to the current content page.




Figure 1- Breadcrumb

Edit to Click Path

To improve navigation the APS menu has changed. The user will select the desired software package i.e., Lab Analytics, PBM, Outreach, Lab Stewardship. Available menu items will then be displayed by selecting the 3 vertical line icon. Administrative settings are located under the settings icon (cog) . Manage data is accessed by selecting the double arrow icon.




Figure 2- Click Path Enhancement

2. Excel File Download

Data mapping download extensions have been updated from .csv to .xls.



Figure 3- Administrative and Manage Data Icon Changes

3. Laboratory Quality and Service Scorecard and Detail Display

For improved differentiation of metric stoplight coding, metric reporting increased from tenth place. Metrics are now reported to the thousandth place.



Figure 5- Laboratory Quality and Service Scorecard View




Figure 6- Laboratory Quality and Service Detail View

4. Laboratory Quality and Service Detail Date Range

An additional date range option has been added to the Quality and Service Detail page for easier month end reporting. The previous month option will select the dates for the last full month. Example: On May 10 the user selects “previous month”, the content page will filter to April 1 to April 30.




Figure 7- Previous month Date Filter

5. Laboratory Orders Query TAT Metric Update

Additional TAT metrics have been added to the Laboratory > Volumes > LIS Orders Query screen. The median TAT (50th percentile) metric remains and the options of 75th and 90th percentile are available using a dropdown menu.

Order to Collect TAT

- In the Order to Collect TAT segments negative time stamps have been removed from the data set for more accurate TAT reporting.




Figure 8- Additional TAT segments

6. Laboratory CPT Detail

Organization of display

- The CPT detail screen will display the three columns (Facilities, Department and CPT Code) in descending order.

Chart title during filtering

- As the CPT Detail screen is filtered the visualization tile will update according to filtered parameters.



Figure 9- Laboratory CPT Detail Chart Title

7. Message Center (Defect)

Message Center was no longer rendering on selection, issue has been resolved.

Date	From	Subject	Status
05/02/23 10:03 pm	Celena Blackwood	PowerPoint Export for Tableau workbook Productivity Scorecard [Celena Blackwood]	
05/02/23 10:02 pm	Celena Blackwood	PHG Export for Tableau workbook LIS Orders Query [Celena Blackwood]	

Figure 10- Message Center

8. Laboratory Quality and Service Detail Filtering (Defect)

For several APS clients, facility and collector type dropdown menus were not populating when accessing the filter funnel. The issue has been located and resolved.

© Accumen Performance Suite

This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen, Inc.