Product Release Notes

Orders Management V7

Release 7.7.2

Aug 2023

Product Detail and Version

Product Name	Version released	Date Released	Comment
Orders Management Portal	7.7.2	Aug 2023	

Revision History

Date	Version	Description	Author
8/11/2023	1	Release Notes	Lynn Lund

Feature release list contained in Version 7.7.2

- 1. User Login Screen View Entered Password
- 2. User Login Error Messaging
- 3. Fields Added to HL7 Template
- 4. MSPQ Limit Date Range for Retirement Questions

Detailed Product Release Description

1. User Login Screen View Entered Password

Previously the user could not read what was typed into the password field on the OM Portal login page. The new enhancement allows the user to select the eye icon to view the typed password and verify the entry before signing in.

Orders Management	Orders Management
TestUser	TestUser
····· v	Typed Password
Sign In	Sign In
Forgot Username Forgot Passw	vord Forgot Username Forgot Password

Figure 1- Using Eye icon to view entered password on Login Screen

2. User Login Error Messaging

Updates have been made to show discrete onscreen error if a user cannot log in due to incorrect information being entered or their user profile is not active. The changes are to provide a better user experience and assist Admin users addressing reported log in issues.

a. When logging into the portal using an incorrect username, the user will see an onscreen message: Incorrect username. Users can at once repeat the login process by entering the correct information or utilizing the *Forgot Username* link.

Orders Manage	ement
Incorrect userna	ame
lundcli	
Password	U
Sign In	
Forgot Username	Forgot Password

Figure 2- Incorrect username attempt onscreen message

b. When logging into the portal using the correct username, but incorrect password, the user will see an onscreen message: Failed password attempt, this is your X attempt since your last successful log in. The 'X' indicates the number of consecutive invalid password attempts for the user. The user will also receive an email to the email address linked to their user profile with the reason for the error. Users can at once repeat the login process by entering the correct information or utilizing the *Forgot Password* link.

© Accumen Technology Solutions Confidential Information Page 4 of 9 This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen Technology Solutions.

Orders Management

Failed passwore your last succes	d attempt, this is your 1st attempt since ssful log in
lundcli1	
Password	ω
Sign In	
Forgot Username	Forgot Password

Figure 3- 1st Failed password attempt onscreen message

- c. Within the Security settings, the client sets the allowed consecutive number of password failures before a user would be temporarily locked out of the OM Portal application and how long the Temporary suspension should last. For the user:
 - i. When logging into the portal using an incorrect password to exceed the threshold set in the security settings, the user will see an onscreen message: Your account is temporarily suspended due to too many unsuccessful attempts. Please retry in Z minutes. The user will also receive an email to the email address linked to their user profile with the reason for the error and the time they are locked out for.

Your account is unsuccessful at	temporarily suspended due to too many tempts. Please retry in 5 minutes
lundcli1	
Password	ω
Sign In	
Forgot Username	Forgot Password

Figure 4- Failed password attempt, Temporary suspension onscreen message

ii. After the locked out time elapses, Users can repeat the login process by entering the correct information or utilizing the *Forgot Password* link.

Admin users can verify the status by checking the user's security state under the User List. Search for the user profile and click the *Admin* link in the record row. The Admin screen will open and provide the time the account will be locked out until.

© Accumen Technology Solutions Confidential Information Page 5 of 9 This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen Technology Solutions.

	lundcli1		٩	All		•		
U	Jsername	Name	E-mail		Security User	Role(s)	Location	
lu	undcli1	test Practi			LockedOut	Practice Admin	10004, Account9	Admin



- d. Within the Security settings, the client sets the timeframe of inactivity that will cause a profile to become suspended.
 - If a user in the 'Suspended' status tries to log into the portal they will see an onscreen message: Your account is suspended, please contact your administrator. The user will also receive an email to the email address linked to their user profile with the reason for the error.

Your account is administrator	s suspended, please contact your
lundly01	
Password	అ
Sign In	
orgot Username	Forgot Password

Orders Management

Figure 6- Suspended user onscreen message

ii. Admin users can verify the status by checking the user's security state under the User List. Clicking the Admin link in the record row will open the Admin screen. Selecting Unsuspend will remove the suspension and allow the user to log in with previously saved password and security questions. If the user cannot remember their information, or had never completed setting up their user profile, also selecting reset will revert the user to a New state, generating a welcome email that will allow the user to set new information.

User List					Export + A
lundly	User- lundly01 is	Active, Suspended		x	
Username	Username	Send Reminder	Create Date	5/18/2020 12:40:50 PM	
lundly01	Password	Password Reset	Password Set Date	1/18/2022 10:37:01 AM	Admin
	Temporary Password	Create Temporary Password	Temporary Password Set Date	7/4/2022 1:22:04 PM	
	Suspended	UnSuspend	Suspended since	7/4/2022 1:22:04 PM	
	User Reset	Reset	Active since	11/13/2020 6:43:10 PM	
			Locked out until		
			Blocked since		
			Close		

Figure 7- Suspended user Admin menu

- e. When a user profile is created or reset, the user needs to use the link in the generated Welcome email to activate their profile and set a password and security questions. Until those steps are completed, the user is in a 'New' status and cannot log into the portal through the normal route. The links in the Welcome email are only valid for the time set in the Security configuration, typically 72 hours.
 - i. If a user in the 'New' status tries to log into the portal they will see an onscreen message: Your account is not activated, please contact your administrator.

Orders Management

Your account is administrator	not ac	tivated, j	please co	ntact your
lundcli1				
Password	v			

Figure 8- New user onscreen message

 ii. Admin users can verify the status by checking the user's security state under the User List. Clicking the *Admin* link in the record row will open the Admin screen.
Selecting Activate will regenerate the new user Welcome email with the directions and activation link for the user to follow.

lundcli			Q,	All	•			
Username	Name	E-mail			Security User State	Role(s)	Location	
lundcli2	test Lund				Active	Client User - Orders and R	10004, 6910, 9009000493	Admin
lundcli1	test Practi				New	Practice Admin	10004, Account9	Admin



3. Fields Added to HL7 Template

The OM-Portal application contains a configuration setting allowing for the hospital or lab customer to designate the fields and structure in the electronic message sent to the Laboratory Information System. With this release configuration the following fields may be added to the HL7 template upon request:

- a. Healthplan type can be added to HL7 segment IN1-15.
- b. Copy to Clinicians name and fax number can be added to the HL7 segment OBR-28.
- c. Guarentor Employer can be added to GT1-51

4. MSPQ Limit Date Range for Retirement Questions

As part of the Medicare Secondary Payor (MSP) Decision tree, the retirement date may need to be entered. Validation is added to the retirement date question to not allow the entry of a future date. If using the calendar, future dates are greyed out. If typing into the date box, the user will not be able to select yes and continue until a valid date is entered.

MSPQ Questionnaire								×	
Is the patient eligible for Medicare due to age? Yes No Is the Patient Employed? Yes No Is the Patient retired? Yes No Enter retirement date and select Yes to continue mm/dd/yyyy Yes Yes									
Reset View MSPQ Save Cancel	Augu	ust 20	23 -			\uparrow	\downarrow		
	Su	Мо	Tu	We	Th	Fr	Sa	~	
Bomovo Cupro	30	31	1	2	3	4	5		
Renove Guara	6	7	8	9	10	11	12		
	13	14	15	16	17	18	19		
Self • Stre	20	21	22	23	24	25	26		
Latest Stre	27	28	29	30	31	1	2		
Lytost	3	4	5	6	7	8	9		
	Cle	ear				То	day		

Figure 10– MSPQ Retirement question calendar with future dates greyed out

© Accumen Technology Solutions Confidential Information Page 8 of 9 This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen Technology Solutions.

User List

© Accumen Technology Solutions Confidential Information Page 9 of 9 This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen Technology Solutions.