

# MyBloodHealth 3.2.2



## Product Release Notes

MyBloodHealth®  
Release 3.2.2  
September 17, 2024

### Product Detail and Version

Product Name	Version	Date	Comment
MyBloodHealth	3.2.2	September 17, 2024	

### Revision History

Date	Version	Description	Author
09.12.2024	1	Release Notes	Tiffany Hall
09.16.2024	1.1	Update - Anemia Treatment Enhancement	Tiffany Hall

### Feature Release list contained in Version 3.2.2

#### New + Improved

1. Change patient height from meter to cm in patient detail section
2. Allow the same user (email address) to be a part of multiple accounts
3. Enhance patient page to display whether the visit was surgical or non-surgical
4. Framework to Show Anemia Treatment on Patient Workspace for Rhapsody Customers
5. Enhance MFA "Verify your identity" instructions to remove the word "one-time" for password application authenticator
6. Add option to delete comment
7. Adding Ferritin level to Worklist

# Detailed Product Release Description

## 1. Change patient height from meter to cm in patient detail section


Previously when a user looked up a patient height from the patient details section it displayed in meters. An enhancement was made display the height in center meters.

## 2. Allow the same user (email address) to be a part of multiple accounts

Previously a user email address could only be assigned to one client account. An enhancement was made to allow an MBH Admin to assign a user to multiple accounts based on the users' job responsibilities by client contract.

## 3. Enhance patient page to display whether the visit was surgical or non-surgical


Previously when a user viewed the patient page it did not indicate if the visit was a surgical or non-surgical visit. An enhancement was made to display the specific visit type of surgical or non-surgical.



## 4. Framework to Show Anemia Treatment on Patient Workspace for Rhapsody Customers


*(Functionality will need to be tested post release in production environment with our Rhapsody customers).*

Previously when a user viewed the patient visit page for a patient who received an anemia management medication therapy it did not display the specific treatment provided for optimization. An enhancement was made to display the specific anemia medication ordered, and the date the treatment was provided to the patient if it was ordered within 90 days from the procedure/visit date.



## 5. Enhance MFA "Verify your identity" instructions to remove the word "one-time" for password application authenticator

Previously when a user accessed MyBloodHealth the multi-factor authenticator instructions were not clear and indicated the password code process was a one-time process. An enhancement was made to remove the word "one-time" from the instructions.



---


## 6. Add option to delete comment

Previously when a user made an error in the comment section there was no functionality to delete the comment. An enhancement was made to add a button to delete a comment.

---


## 7. Adding Ferritin level to Worklist

Previously when a user viewed the worklist of patients it did not display a ferritin result if the patient had a result reported. An enhancement was made to display the ferritin result if the test was reported within 60 days from the patient procedure/visit date. The result value will display in black font if the result is normal and if the result is abnormal it will display in red font.



---

## Downloads



MBH-MyBlood... 226.pdf  
13 Sep 2024, 11:33 AM

---

© MyBloodHealth, an Accumen Company

This document represents proprietary information and neither it nor its contents may be disclosed to anyone without the specific permission of Accumen, Inc.